

Lumbar Puncture Trainer Set Up Procedure


Antonio Sessoms - 2024-08-28 - Low-Fidelity

Description:

Lumbar puncture and epidural anesthesia trainer. Intended for single or team training scenarios. Includes five types of puncture blocks to enhance training with different types of challenges.

Equipment:

1 lumbar region model / 5 lumbar puncture blocks types : 2 normal CSF, 1 obesity CSF, 1 senior CSF, 1 senior obesity CSF, 1 epidural / 1 lumbar region skin cover / 1 supporting stand/ 1 support base for team teaching / 1 lumbar spine model / 1 irrigation bag / 1 syringe / 1 stand for irrigation bag / 1 instruction manual / 1 storage case



| | | | |
|--|---|---------------------------------------|---|
| a. Lumbar Region Model | 1 | e. Skin Cover | 1 |
| b. CSF Puncture Blocks (4 variations 5 pieces) | | f. Syringe | 1 |
| b-1. Normal CSF | 2 | g. Irrigator Pole | 1 |
| b-2. Obesity CSF | 1 | h. Irrigator Bag | 1 |
| b-3. Senior CSF | 1 | i. Supporting Stand (Two Way) | 1 |
| b-4. Senior Obesity OSF | 1 | j. Support Base (Team Teaching) | 1 |
| c. Epidural Puncture Block | 1 | | |
| d. Lumbar Spine Model | 1 | | |

Set Up Instructions:

1. Position the Supporting Stand
2. Hang irrigator bag on irrigator pole
3. Select irrigation block
 1. "Note arrow marking on base of puncture block"
4. Connect puncture block to irrigator bag (on end with arrow)
5. Connect a syringe to opposite end of block. Insert and turn clockwise and lock in syringe in place
6. Release the clamp
7. Fill puncture block with water. Tilt the puncture block and aspirate water with syringe until only water remains in the tubing
8. Close the clamp
9. Attach back skin to Lumbar Region Model, noting the marks on inner skin cover with arrow pointing towards the tubing of the irrigator bag